


Antigrade Femoral Nail

(Short & Long)
(S.S. & Titanium)

SMPL™



Proximal diameter 12.5mm

Cannulated self-tapping 6.4mm
Leg screw length 70 to 120mm (5mm Variation)

Cannulated self-tapping 6.4mm
Leg screw length 70 to 120mm (5mm Variation)

CCD Angle 135°

Anatomical Angle 6°degree

Distal diameter of 9, 10, 11 & 12mm

Locking Screw 4.9mm
Length 25 to 70mm

Dynamization 5mm



◀◀ SMPL Antigrade Femoral Nail

Code No. S182 / S182T (Dia. 9, 10, 11 & 12)
Stainless Steel & Titanium

SHORT NAIL

Code No. S182-Dia-240 (Stainless Steel)
Code No. S182T-Dia-240 (Titanium)

LONG NAIL (L & R)

Code No. S182-Dia(L&R)-340 To 440 (20mm Variation) (Stainless Steel)
Code No. S182T-Dia(L&R)-340 To 440 (20mm Variation) (Titanium)



Leg Screw Ø 6.4 mm (Cannulated) ▶▶▶

Stainless Steel Code No. S175-70 To 120 mm (5mm Variation)
Titanium Code No. S175T-70 To 120 mm (5mm Variation)



Locking Screw Ø 4.9 mm ▶▶▶

Stainless Steel
Code No. S155-25 To 70 mm (5mm Variation)
Titanium
Code No. S155T-25 To 70 mm (5mm Variation)

End Cap ▶▶▶

Stainless Steel
Code No. S183-(2,5,10,15 & 20)
Titanium
Code No. S183T-(2,5,10,15 & 20)



All. Sighting Zig-proximal & Short Distal ▶▶▶

Cat. No. : S255-01



Nail Insertion Handle ▶▶▶

Cat No. S255-03



Top Nut for Sighting Zig ▶▶▶

Cat No. S255-04



Conical Bolt ▶▶▶

Cat No. S255-05



Box Spanner for Conical Bolt ▶▶▶

Cat No. S255-06



Stop Drill Sleeve for Proximal Reamer ▶▶▶

Cat No. S255-07



Reamer Dia 13.0mm Cannulated ▶▶▶

Cat No. S255-08



Reamer Dia 6.4mm Cannulated of Leg Screw ▶▶▶

Cat No. S255-09



Cylindrical Cross Head Bolt ▶▶▶

Cat No. S255-11



Socket Head Bolt ▶▶▶

Cat No. S255-12



T-Handle for Socket Head Bolt ▶▶▶

Cat No. S255-13



Sleeve for Proximal Leg Screw ▶▶▶

Cat No. S255-14



Protection Sleeve ▶▶▶

Cat. No. : S255-15



Guide Wire Sleeve Dia 2.2mm ▶▶▶

Cat No. S255-16



Drill Sleeve Dia 4mm ▶▶▶

Cat No. S255-17



Depth Gauge ▶▶▶

Cat No. S255-18



Extraction Rod with Sliding Hammer ▶▶▶

Cat No. S255-20



Hammer ▶▶▶

Cat No. S255-21



Cannulated Screw Driver for Leg Screw 6.4mm ▶▶▶

Cat No. S255-22



Locking Screw Driver 4.9mm ▶▶▶

Cat No. S255-23



Chuck ▶▶▶

Cat No. S255-24



Drill Bit 4mm ▶▶▶

Cat No. S255-25



Guide Pin ▶▶▶

Cat No. S255-26



Trocar ▶▶▶

Cat No. S255-27



SMIT MEDIMED PVT. LTD.

Plot No.10, Phase-1, B/h. Prashant Eng., G.I.D.C. Vatva, Ahmedabad-382 445. (Gujarat) INDIA.

Phone: +91-79-25896751, 25896752 // Fax : +91-79-25894150

E-mail : smitmed1@sancharnet.in // Visit us : www.smitmedimed.com